

The Cor.At Starter Kit enables the customer to perform a quick and easy initial experiment in his field of research. The kit includes a 500k vial of Cor.At GFP cardiomyocytes, Cor.At culture media as well as all necessary disposables needed for thawing and seeding of the cells.

The flexible design of the kit allows for various applications, e.g. electrophysiology, toxicity screening, and immunostaining. The Cor.At starter kit helps the new customer to feel comfortable in handling Cor.At cardiomyocytes before starting key experiments.

Cor.At Flowchart

